

K1826 蛋白抗体

产品货号： mlR17215

英文名称： K1826

中文名称： K1826 蛋白抗体

别名： MSANTD4; Coiled coil domain containing protein KIAA1826; DC25; Hypothetical protein LOC84437; KIAA1826; Myb/SANT-like DNA-binding domain containing 4 with coiled-coils; Myb/SANT-like DNA-binding domain-containing protein 4.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Rabbit, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需

做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 41kDa

细胞定位 : 细胞核

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human K1826:1-100/345

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : The function of this protein has yet to be defined.

Similarity:

Contains 1 Myb-like domain.

SWISS:

Q8NCY6

Gene ID:

84437

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.