

整联蛋白重复蛋白 2 抗体

产品货号： mlR18752

英文名称： ITFG2/MDS028

中文名称： 整联蛋白重复蛋白 2 抗体

别名： Integrin alpha FG GAP repeat containing 2; Integrin alpha FG GAP repeat containing protein 2; Integrin-alpha FG-GAP repeat-containing protein 2; ITFG 2; ITFG2; ITFG2_HUMAN; MDS 028; MDS028; OTTHUMP00000161190; RGD1311944; Uncharacterized hematopoietic stem/progenitor cells protein MDS028.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Chicken, Dog, Pig, Horse, Rabbit, Sheep, Guinea Pig, Chimpanzee, Rhesus monkey, Gorilla, Orangutan

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 49kDa

细胞定位： 细胞核 细胞浆

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human ITFG2/MDS028:201-300/447

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : ITFG2 (Integrin-alpha FG-GAP repeat-containing protein 2) is a 447 amino acid protein that contains two FG-GAP repeats, a motif commonly found in integrin proteins.

Similarity:

Contains 2 FG-GAP repeats.

SWISS:

Q969R8

Gene ID:

55846

Important Note:



This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.