

锚蛋白重复结构域蛋白 32 抗体

产品货号： mlR9749

英文名称： ANKRD32

中文名称： 锚蛋白重复结构域蛋白 32 抗体

别名： Brctx; Ankrd32; Ankyrin repeat domain-containing protein 32; ANR32_HUMAN; BRCT domain-containing protein 1; Brctd1.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Horse, Rabbit,

产品应用： ELISA=1:500-1000 IHC-F=1:400-800 IF=1:50-200 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 121kDa

细胞定位 : 细胞核 细胞浆

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human ANKRD32:21-120/1058

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 Similarity :

Contains 3 ANK repeats.

Contains 2 BRCT domains.

SWISS:

Q9BQI6

Gene ID:

84250

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.