

纳曲酮抗体 IgG

产品货号： mlR0903

英文名称： Natrexone

中文名称： 纳曲酮抗体 IgG

产品类型： 药物与化合物抗体

研究领域： 药物及化合物

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Natrexone

产品应用： ELISA=1:500-1000

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated Natrexone:

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.



保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed

产品介绍： Naltrexone is a drug prescribed to assist people in maintaining abstinence after they have detoxified or withdrawn from opioids (eg heroin) or alcohol. Naltrexone is classed as an opioid agonist. It works by blocking the opioid receptors in the brain so that opiates have no effect on the brain.

SWISS:

N/A

CAS:

16590-41-3

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.