

## 丝氨酸/苏氨酸蛋白激酶 32A 抗体

产品货号： mIR12832

英文名称： STK32A

中文名称： 丝氨酸/苏氨酸蛋白激酶 32A 抗体

别 名： A930015B13Rik; MGC22688; Serine/threonine kinase 32A; Serine/threonine protein kinase 32A; STK 32A; YANK 1; YANK1; Yet another novel kinase 1.

研究领域： 细胞生物 信号转导 激酶和磷酸酶

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Cow, Horse, Rabbit, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 46kDa

细胞定位： 细胞膜

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

**免 疫 原：** KLH conjugated synthetic peptide derived from human STK32A:161-260/396

**亚 型：** IgG

**纯化方法：** affinity purified by Protein A

**储 存 液：** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed

**产品介绍：** STK32A is a member of the serine/threonine protein kinase family. The specific function of this kinase is not known.

**Subcellular Location:**

Cell membrane.

**Similarity:**

Belongs to the protein kinase superfamily.

Ser/Thr protein kinase family.

Contains 1 protein kinase domain.

**SWISS:**

Q8WU08

**Gene ID:**

202374

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.