

HRASLS5 蛋白抗体

产品货号： mIR17386

英文名称： HRASLS5

中文名称： HRASLS5 蛋白抗体

别名： 1700055A14Rik; 1700086E08Rik; 4921526K24Rik; H-rev107-like protein 5; HRAS like suppressor family, member 5; HRAS-like suppressor 5; HRASLS5; HRSL5_HUMAN; iNAT; MGC116140; RGD1308376; RLP1.

研究领域： 肿瘤 细胞生物 细胞类型标志物

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 30kDa

细胞定位 : 细胞浆

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human HRASLS5:61-160/279

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.



PubMed : PubMed