

## 庆大霉素抗体

产品货号： mlR0972

英文名称： Gentamicin

中文名称： 庆大霉素抗体

别 名： Garamycin; Gentamicin sulfate salt; Gentamycin ;Gentiomycin C.

产品类型： 药物与化合物抗体

研究领域： 药物及化合物

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Gentamicin

产品应用： ELISA=1:500-1000

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated to Gentamicin:

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed

**产品介绍：** Gentamicin is an aminoglycoside antibiotic, used to treat many types of bacterial infections, particularly those caused by Gram-negative bacteria. However, gentamicin is not used for Neisseria gonorrhoeae, Neisseria meningitidis or Legionella pneumophila bacterial infections (because of the risk of the patient going into shock from lipid A endotoxin found in certain gram negative organisms).

**SWISS:**

N/A

**CAS:**

1403-66-3

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.