

肿瘤睾丸抗原 1B1A

产品货号： mlR4501

英文名称： testis antigen 1B

中文名称： 肿瘤/睾丸抗原 1B/1A

别名： Autoimmunogenic cancer/testis antigen NY ESO 1; testis antigen 1A; Autoimmunogenic cancer/testis antigen NY-ESO-1; Cancer antigen 3; Cancer/testis antigen 1; Cancer/testis antigen 1B; Cancer/testis antigen 6.1; CT6.1; CTAG 1; CTAG 1B; CTAG; CTAG1; CTAG1A; CTAG1B; CTG1B_HUMAN; ESO 1; ESO1; L antigen family member 2; LAGE 2; LAGE 2 protein; LAGE 2B; LAGE-2; LAGE2; LAGE2 protein; LAGE2A; LAGE2B; New York esophageal squamous cell carcinoma 1; NY ESO 1; NYESO 1; NYESO1.

研究领域： 肿瘤 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:50-200 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 18kDa

细胞定位 : 细胞浆

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human Cancer/testis antigen 1B/1A:75-180/180

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed : PubMed

产品介绍： Expressed in testis and ovary and in a wide variety of cancers. Detected in uterine myometrium. Expressed from 18 weeks until birth in human fetal testis. In the adult testis, is strongly expressed in spermatogonia and in primary spermatocytes, but not in post-meiotic cells or in testicular somatic cells (at protein level).

Subcellular Location:

Cytoplasm.

Tissue Specificity:

Expressed in testis and ovary and in a wide variety of cancers. Detected in uterine myometrium. Expressed from 18 weeks until birth in human fetal testis. In the adult testis, is strongly expressed in spermatogonia and in primary spermatocytes, but not in post-meiotic cells or in testicular somatic cells (at protein level).

Similarity:

Belongs to the CTAG family.

SWISS:

P78358

Gene ID:

1485

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

