

人类疱疹病毒 8 抗体

产品货号： mlR0746

英文名称： HHV8 ORF50

中文名称： 人类疱疹病毒 8 抗体

别 名： HHV 8; Human Herpes Virus 8; Human herpesvirus 8; HHV8 ORF50; Kaposi's sarcoma associated herpes virus; Kaposi's sarcoma associated herpesvirus; KSHV; ORF50 protein.

研究领域： 肿瘤 免疫学 细菌及病毒

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： HHV8

产品应用： ELISA=1:500-1000 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

分 子 量： 74kDa

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from HHV8 ORF50:2-100/691

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed

产品介绍： HHV8 has been found to be associated with three different diseases observed in AIDS patients; kaposi's sarcoma, primary effusion lymphoma (which is a rare type of non-Hodgkin lymphoma affecting the body cavities) and multicentric Castleman's disease. To date there is much evidence to support a direct role for HHV8 in kaposi's sarcoma.

SWISS:

N/A

Gene ID:

4961526

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.