

羧基端小结构域磷酸酶 2 抗体

产品货号： mlR14099

英文名称： CTDSP2

中文名称： 羧基端小结构域磷酸酶 2 抗体

别名： Carboxy terminal domain RNA polymerase II polypeptide A small phosphatase 2; Carboxy-terminal domain RNA polymerase II polypeptide A small phosphatase 2; CTDS2_HUMAN; Ctdsp2; NIF2; NLI interacting factor 2; NLI-interacting factor 2; Nuclear LIM interactor interacting factor 2; Nuclear LIM interactor-interacting factor 2; OS4; Protein OS 4; Protein OS-4; PSR2; SCP2; Small C-terminal domain phosphatase 2; Small CTD phosphatase 2.

研究领域： 肿瘤 细胞生物 染色质和核信号 信号转导 肿瘤细胞生物标志物 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig, Cow, Horse, Rabbit, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 31kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human CTDSP2:25-150/271

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed： PubMed

产品图片

