

## F-box 蛋白相关蛋白 6 抗体

产品货号： mlR6415

英文名称： FbxO6

中文名称： F-box 蛋白相关蛋白 6 抗体

别名： F box protein that recognizes sugar chains 2; F box/G domain protein 2; F box G domain protein 2; F box only protein 6; F box protein 6; F box protein FBG2; F box protein Fbx6; FBG 2; FBG2; FBS 2; FBS2; FBX 6; FBX6; FBXO 6.

研究领域： 细胞生物 免疫学 信号转导

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Cow,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量：34kDa

细胞定位：细胞浆

性状：Lyophilized or Liquid

浓度：1mg/ml

免疫原：KLH conjugated synthetic peptide derived from human FbxO6:101-200/293

亚型：IgG

纯化方法：affinity purified by Protein A

储存液：0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件：Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed：PubMed

产品介绍：The F box, named after cyclin F in which it was originally observed, is an approximately 40-amino acid motif that binds SKP1. F-box proteins are components of modular E3 ubiquitin protein ligases called SCFs

(SKP1, cullin, F-box proteins), which function in phosphorylation-dependent ubiquitination.

**Function:**

Substrate-recognition component of some SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complexes. Involved in endoplasmic reticulum-associated degradation pathway (ERAD) for misfolded luminal proteins by recognizing and binding sugar chains on unfolded glycoproteins that are retrotranslocated into the cytosol and promoting their ubiquitination and subsequent degradation. Able to recognize and bind denatured glycoproteins, which are modified with not only high-mannose but also complex-type oligosaccharides. Also recognizes sulfated glycans. Also involved in DNA damage response by specifically recognizing activated CHEK1 (phosphorylated on 'Ser-345'), promoting its ubiquitination and degradation. Ubiquitination of CHEK1 is required to insure that activated CHEK1 does not accumulate as cells progress through S phase, or when replication forks encounter transient impediments during normal DNA replication.

**Subunit:**

Interacts with VCP.

**Subcellular Location:**

Cytoplasm.

**Similarity:**

Contains 1 F-box domain.

Contains 1 FBA (F-box associated) domain.

**SWISS:**

Q9NRD1

**Gene ID:**

26270

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

