

肥大细胞蛋白酶 4 抗体

产品货号： mlR18733

英文名称： MCPT4

中文名称： 肥大细胞蛋白酶 4 抗体

别 名： Mast cell protease 4; Mcp 4; Mcpt4; MCPT4_MOUSE; MMCP 4; MMCP 4A; MMCP 4B; MMCP-4; MMCP4; MMCP4A; MMCP4B; MSMCP; Myonase; Serosal mast cell protease.

研究领域： 细胞生物 免疫学 泛素

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Mouse, Rat,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500

（石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 25kDa

细胞定位： 细胞浆 细胞膜

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from mouse MCPT4:51-150/246

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed