

## m Taq 酶抗体

产品货号： mIR10636

英文名称： m Taq

中文名称： m Taq 酶抗体

别 名： Taq DNA Polymerase

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Taq

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500  
(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 91kDa

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : Taq DNA Polymerase:

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20 ° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed :** PubMed

#### 产品介绍 **Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.