

LYRM1 蛋白抗体

产品货号： mlR18582

英文名称： LYRM1

中文名称： LYRM1 蛋白抗体

别名： 1110065L10Rik; 2310004B22Rik; 4930404J24Rik; A211C6.1; LYR motif containing 1; LYR motif containing protein 1; LYR motif-containing protein 1; LYRM 1; lyrm1; LYRM1_HUMAN; OTTHUMP00000162292; OTTMUSP00000031126; OTTMUSP00000031127; OTTMUSP00000031128; OTTMUSP00000031129; RGD1563498.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Cow, Horse, Rabbit, Sheep,

产品应用： IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 14kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human LYRM1:1-80/122

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20 ° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed