

## 磷酸化肌肉型 6 磷酸果糖激酶/磷酸果糖激酶 1 抗体

产品货号： mlR16457

英文名称： phospho-PFKM/PFK1 (Ser775)

中文名称： 磷酸化肌肉型 6 磷酸果糖激酶/磷酸果糖激酶 1 抗体

别名： PFKM (phospho S775); PFKM (phospho Ser775); p-PFKM (Ser775); Fructose 6 Phosphate Kinase; 6 Phosphofructokinase Muscle Type; GSD7; PFKA; PFK-A; PFKP; PFKX; Phosphofructo 1 Kinase Isozyme A; Phosphofructo-1-kinase isozyme A; Phosphofructokinase 1; Phosphofructokinase M; Phosphofructokinase-M; Phosphofructokinase, muscle; Phosphofructokinase, muscle type; Phosphofructokinase, polypeptide X; Phosphohexokinase; PFKAM\_HUMAN; 6-phosphofructokinase, muscle type; PFK1; PFK-1; ATP-PFK; PPP1R122.

产品类型： 磷酸化抗体

研究领域

抗体来源： Rabbit

克隆类型： Polyclonal



交叉反应： Human, Mouse, Rat, Pig, Cow, Horse, Sheep,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500  
(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 86kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated Synthesised phosphopeptide derived from human PFKM around the phosphorylation site of Ser775:KR(p-S)GE

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed

**产品介绍 background:**

Three phosphofructokinase isozymes exist in humans: muscle, liver and platelet. These isozymes function as subunits of the mammalian tetramer phosphofructokinase, which catalyzes the phosphorylation of fructose-6-phosphate to fructose-1,6-bisphosphate. Tetramer composition varies depending on tissue type. This gene encodes the muscle-type isozyme. Mutations in this gene have been associated with glycogen storage disease type VII, also known as Tarui disease. Alternatively spliced transcript variants have been described.[provided by RefSeq, Nov 2009]

**Function:**

Catalyzes the third step of glycolysis, the phosphorylation of fructose-6-phosphate (F6P) by ATP to generate fructose-1,6-bisphosphate (FBP) and ADP.

**Subunit:**

Homo- and heterotetramers. Muscle is M4, liver is L4, and red cell is M3L, M2L2, or ML3. Interacts (via C-terminus) with HK1 (via N-terminal spermatogenic cell-specific region)

**Subcellular Location:**

Cytoplasm.

**Post-translational modifications:**

GlcNAcylation decreases enzyme activity.

**DISEASE:**

Glycogen storage disease 7 (GSD7) [MIM:232800]: A metabolic disorder characterized by exercise intolerance with associated nausea and vomiting, muscle cramping, exertional myopathy and compensated hemolysis. Short bursts of intense activity are particularly difficult. Severe muscle cramps and myoglobinuria develop after vigorous exercise.

**Similarity:**

Belongs to the phosphofructokinase type A (PFKA) family. ATP-dependent PFK group I subfamily. Eukaryotic two domain clade 'E' sub-subfamily.

**SWISS:**

P08237

**Gene ID:**

5213

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

