

KTI12 蛋白抗体

产品货号： mlR16862

英文名称： KTI12

中文名称： KTI12 蛋白抗体

别名： KTI12 homolog chromatin associated (S. cerevisiae); KTI12 homolog chromatin associated; MGC20419; Protein KTI12 homolog; RP11 91A18.3; SBB181; TOT4. KTI12_HUMAN

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Horse,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 39kDa

细胞定位 : 细胞外基质

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human KTI12:261-354/354

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : background:

The function of KTI12 is not known.

Similarity:

Belongs to the KTI12 family.

SWISS:

Q96EK9

Gene ID:

112970

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.