

蛋白酪氨酸磷酸酶样 A 结构域蛋白 2 抗体

产品货号： mIR19588

英文名称： PTPLAD2

中文名称： 蛋白酪氨酸磷酸酶样 A 结构域蛋白 2 抗体

别名： Protein tyrosine phosphatase like A domain containing protein 2.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500

（石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 28kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human PTPLAD2:1-100/232 <Cytoplasmic>

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : PTPLAD2 is a protein that belongs to the PTPLA family. The function of PTPLAD2 remains unknown.

Subcellular Location:

Membrane; Multi pass membrane protein

SWISS:

Q5VWC8

Gene ID:

401494

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.