

## 小牛胸腺 DNA 抗体

产品货号： mlR0313

英文名称： calf thymus DNA

中文名称： 小牛胸腺 DNA 抗体

别名： DNA,Calf ThyMus,Na; dna single stranded;DNA from calf thymus;DEOXYRIBONUCLEIC ACI;  
CalfthymusDNA pyrolyzate; DNABOVINE (CALF) THYMUS; DNA, Na Salt, Calf Thymus; Deoxyribonucleic acids,  
thymus; DEOXYRIBONUCLEIC ACID, GENOMIC; Deoxyribonucleicacid(calfthymus)

研究领域： 细胞凋亡

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Cow,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： Deoxyribonucleic acid sodium salt from calf thymus:full length

亚型： IgG

纯化方法： affinity purified by Protein A



**储 存 液：** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed