

## 乙型肝炎病毒 X 蛋白 HBX 抗体

产品货号： mlR2146

英文名称： Hepatitis B virus pre-X protein

中文名称： 乙型肝炎病毒 X 蛋白 HBX 抗体

别 名： HBX; pre-X protein; HB-X; pre X protein; HBV X protein; pre-X protein [Hepatitis B virus].

研究领域： 细菌及病毒

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Hepatitis B virus

产品应用： ELISA=1:500-1000

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 23kDa

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from Hepatitis B virus X protein:1-50/210

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed

**产品介绍：** Hepatitis B virus

**SWISS:**

N/A

**Gene ID:**

N/A

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.