

跨膜蛋白 9 超家族成员 3 抗体

产品货号： mIR19944

英文名称： TM9SF3

中文名称： 跨膜蛋白 9 超家族成员 3 抗体

别名： Endomembrane protein emp70 precursor isolog; EP70 P iso; EP70-P-iso; RGD1564625; RP1134E5.1; SM 11044 binding protein; SM-11044-binding protein; SMBP; TM9S3_HUMAN; Tm9sf3; Transmembrane 9 superfamily member 3; Transmembrane 9 superfamily protein member 3; UNQ245/PRO282.

研究领域： 细胞生物 免疫学 跨膜蛋白

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Chicken, Dog, Cow, Horse, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 65kDa

细胞定位： 细胞膜

性状： Lyophilized or Liquid

浓度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human TM9SF3:31-130/589

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed