

跨膜蛋白 TMBIM1 抗体

产品货号： mlR19947

英文名称： TMBIM1

中文名称： 跨膜蛋白 TMBIM1 抗体

别名： LFG3; MST100; MSTP100; PP1201; Protein lifeguard 3; Protein RECS1 homolog; RECS1; TMBI1_HUMAN; Tmbim1; Transmembrane BAX inhibitor motif containing 1; Transmembrane BAX inhibitor motif-containing protein 1.

研究领域： 细胞生物 细胞膜受体 结合蛋白 跨膜蛋白 细胞膜蛋白

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Cow, Horse,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 35kDa

细胞定位： 细胞膜

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human TMBIM1:221-311/311

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed