

蚓激酶抗体

产品货号： mlR0484

英文名称： lumbrokinase

中文名称： 蚓激酶抗体

别名： lumbrokinase

研究领域： 激酶和磷酸酶

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Eisenia fetida

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 26kDa

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from Eisenia fetida of lumbrokinase:1-100/245

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.



保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : No data available.

SWISS:

N/A

Gene ID:

N/A

Database links:

GenBank: ABA43718.1

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.