

唾液酸转移酶 8D 抗体

产品货号： mlR17479

英文名称： SIA8D

中文名称： 唾液酸转移酶 8D 抗体

别名： 8-sialyltransferase 8D; 8-sialyltransferase; Alpha 2 8 sialyltransferase 8D; Alpha-2; CMP N acetylneuraminate poly alpha 2 8 sialyl transferase; CMP-N-acetylneuraminate-poly-alpha-2; Polysialyltransferase 1; Polysialyltransferase-1; PST; PST1; SIA8D_HUMAN; Sialyltransferase 8 alpha 2 8 polysialyltransferase D; Sialyltransferase 8D; Sialyltransferase St8Sia IV; SIAT8-D; SIAT8D; ST8 alpha N acetyl neuraminide alpha 2 8 sialyltransferase 4; ST8Sia IV; ST8SIA4; ST8SiaIV.

研究领域： 细胞生物 细胞类型标志物

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Chicken, Cow, Horse, Rabbit, Guinea Pig, African Green Monkey

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 41kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human SIA8D:101-200/359

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.



PubMed : PubMed