

防御素 α 6 抗体

产品货号： mlR7129

英文名称： Defensin alpha 6

中文名称： 防御素 α 6 抗体

别名： alpha 6; alpha 6; DEF6; DEF6_HUMAN; DEFA 6; DEFA6; Defensin 6; Defensin alpha 6; Defensin alpha 6 Paneth cell specific; Defensin; Defensin-6; HD 6; HD6.

研究领域： 细胞生物 免疫学 细菌及病毒

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Mouse, Rat,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 11kDa

细胞定位： 细胞浆 细胞外基质

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from mouse Defensin alpha 6:41-140/100

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed : PubMed

产品介绍 : Has very low antimicrobial activity against Gram-negative and Gram-positive bacteria. May protect cells against infection with HIV-1.

Function:

Has very low antimicrobial activity against Gram-negative and Gram-positive bacteria. May protect cells against infection with HIV-1.

Subunit:

Homodimer.

Subcellular Location:

Secreted.

Tissue Specificity:

Paneth cells of the small intestine.

Similarity:

Belongs to the alpha-defensin family.

SWISS:

Q01524

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

