

高尔基体自身蛋白 8B 抗体

产品货号： mlR16257

英文名称： GOLGA8B

中文名称： 高尔基体自身蛋白 8B 抗体

别名： GOLGA5; Golgi autoantigen golgin subfamily a 8B; GOG8B_HUMAN; Golgin 67; Golgin A8 family member B; Golgin subfamily A member 8B.

研究领域： 细胞生物 信号转导 转运蛋白

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500

(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 67kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human GOLGA8B:511-603/603

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 癆 for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20 癆. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 癆.

PubMed : PubMed

产品介绍 : Function: May be involved in maintaining Golgi structure.

Subcellular Location:

Golgi apparatus ?Golgi stack membrane; Peripheral membrane protein

SWISS:

A8MQT2

Gene ID:

440270

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

