

黑色素瘤相关抗原 H1 抗体

产品货号： mlR18628

英文名称： MAGEH1

中文名称： 黑色素瘤相关抗原 H1 抗体

别名： Apoptosis related protein 1; Apoptosis Related Protein; Apoptosis-related protein 1; APR 1; APR-1; Apr1; MAGE H1 antigen; MAGE-H1 antigen; MAGEH; MAGEH1; MAGH1_HUMAN; melanoma antigen, family H, 1 protein; Melanoma associated antigen H1; Melanoma-associated antigen H1; Restin.

研究领域： 肿瘤 细胞生物

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Horse, Rabbit, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 24kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原 : KLH conjugated synthetic peptide derived from human MAGEH1:31-130/219

亚型 : IgG

纯化方法 : affinity purified by Protein A

储存液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : MAGE-H1 is a 219 amino acid protein that contains a type II MAGE homology domain (MHD). Enhanced ligand stimulation promotes MAGE-H1 interaction with the type II death domain of NGFR p75. It is suggested that MAGE-H1 accelerates differentiation in response to nerve growth factor in cells.

Similarity:

Contains 1 MAGE domain.

SWISS:

Q9H213

Gene ID:

28986

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic



applications.