

兔抗 UBE2U 多克隆抗体

中文名称: 兔抗 UBE2U 多克隆抗体
英文名称: Anti-UBE2U rabbit polyclonal antibody
相关类别: 一抗
抗 原: UBE2U
储 存: 冷冻 (-20℃)
宿 主: Rabbit
反应种属: Human
标 记 物: Unconjugate
克隆类型: rabbit polyclonal

包装:

25 ul 400
100ul 1500
200 ul 2500

技术规格

Background:	UBE2U (ubiquitin-conjugating enzyme E2U) is a 321 amino acid E2 ubiquitin conjugating enzyme that catalyzes the covalent attachment of ubiquitin to other proteins. Existing as two alternatively spliced isoforms, UBE2U is encoded by a gene located on human chromosome 1, which spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome.
Applications:	ELISA, IHC
Name of ant	UBE2U

ibody:	
Immunogen:	Fusion protein of human UBE2U
Full name:	ubiquitin-conjugating enzyme E2U (putative)
SwissProt:	Q5VVX9
ELISA Recom mended dilu tion:	5000-10000
IHC positive control:	Human breast cancer
IHC Recomm end dilution:	25-100