

---

# Recombinant Human IL-10

Catalog Number : mly02501

## Characteristics:

Source	Chinese Hamster Ovary cell line
Description	Human IL-10(Ser19-Asn178) Accession #P22301
Predicted molecular mass	18.6kDa

## Specification:


Appearance	White powder,colorless clear liquid after reconstitution
Purity	≥95%, by SDS-PAGE (under reducing (R)& Non-reducing conditions, visualized by Coomassie staining)
Endotoxin	≤10 EU/mg by the LAL method
Activity	Measured in a cell proliferation assay using MC-9 mouse mast cells. The ED <sub>50</sub> is ≤1ng/ml.
Formulation	Lyophilized from a 0.22 μm-filtered solution containing PBS,5% Mannitol,0.01% Tween 80 , pH7.4
Size	10 μg/vial, 50 μg/vial, 1 mg/vial

## Handling and Storage:

Reconstitution	It is recommended to redissolve in sterile deionized water.
Shipping	The product is shipped with blue ice
Storage & Stability	36 months at -20°C to -80°C in lyophilized state 6 months at -20°C to -80°C under sterile conditions after reconstitution 7-10 days at 2°C to 8°C under sterile conditions after reconstitution Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

## Data:

### Bioactivity



Measured in a cell proliferation assay using MC-9 mouse mast cells.