
Recombinant Human IGF-II (HEK293-expressed)

Catalog Number : mly12441H

Characteristics:

Source	HEK293
Description	Human IGF-II(Ala25-Glu91) Accession # P01344.1
Predicted molecular mass	7.5kDa

Specification:

Appearance	White powder,Colorless clear liquid after reconstitution
Purity	≥95%, by SDS-PAGE (under conditions, visualized by Coomassie staining)
Endotoxin	≤10 EU/mg by the LAL method
Activity	Measured in a serum-free cell proliferation assay using MCF-7 human breast cancer cells. The ED ₅₀ for this effect is ≤50 ng/mL.
Formulation	Lyophilized from a 0.22 μm filtered solution containing TFA,ACN,5% Trehalose.

Handling and Storage:

Reconstitution	It is recommended to redissolve in 4mM HCl.
Shipping	The product is shipped with blue ice
Storage & Stability	36 months at -20°C to -80°C in lyophilized state 6 months at -20°C to -80°C under sterile conditions after reconstitution 7-10 days at 2°C to 8°C under sterile conditions after reconstitution Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Data:

Bioactivity

